



Beryllium Measurement in Commercially Available Wet Wipes

Linda D. Youmans-McDonald, Mike Brisson, Maureen Bernard

Savannah River Nuclear Solutions LLC

Aiken, SC 29808, USA

Anoop Agrawal, John Cronin

Berylliant, Inc.

Tuscon, Arizona

SRNS-N3300-2010-00023



Savannah River

Nuclear Solutions, LLC

A Fluor Daniel Partnership



Disclaimers

- **Mention of commercial products in this presentation does not imply endorsement by the author, SRNS, SRS, or the U.S. Department of Energy (DOE).**
- **The speaker is an employee of a Department of Energy (DOE) contractor, and is not a spokesperson for DOE itself.**





History

- Beryllium sample increase is expected from work generated by American Resource and Recovery Act
- Current Analysis (HF acid digestion and Column Separation, measurement by sequential ICP-ES) is very labor intensive, resulting in long turnaround times.
- No funds are available for additional laboratory personnel or more expensive equipment
- Need emerged for an inexpensive and less labor intensive beryllium measuring method
- Optical Fluorescence Method is less labor intensive
- Costs are significantly less than to purchase an ICP-MS





Challenges

- **Savannah River Site field IH prefer individually packaged pre-wetted wipes to perform beryllium sampling**
- **The ASTM lead wipe standard requires that sample media be wet and individually packaged. A similar beryllium wipe standard is expected to be established in the near future.**
- **Most customers currently using the optical fluorescence method use dry MCE filters**
- **More ammonium bifluoride is required to submerge the commercially available wipes into the solution (detection limit challenges)**





Palintest™, Ghostwipe™, and MCE Blank Testing

	Palintest™ wipe	Ghostwipe™	High Purity™ MCE filter
Blank: 0 µg Be	0.06µg	<0.005µg	<0.005µg
	0.03µg	<0.005µg	<0.005µg
	0.04µg		<0.005µg
			<0.005µg
Spiked with Be acetate			
	0.061µg (153%)	0.04µg (100%)	0.0038µg (93%)
Wipes: 0.04µg Be	0.057µg (139%)	0.041µg (103%)	0.0038µg (93%)
*Filters: 0.004µg	0.064µg (160%)	0.04µg (100%)	0.0038µg (93%)
	0.049µg (120%)	0.037µg (93%)	0.004µg (100%)
	0.063µg (154%)	0.038µg (93%)	0.004µg (100%)

* Filters were spiked with a concentration of beryllium that is less than the established method reporting limit.





Method Reporting Limit for MCE Filters

Standard Conc. Used =	0.08	0.08	0.08
Units	ppb	ppb	ppb
Other Parameters/Information:	filter	filter	filter
Replicate	Measured	Measured	Measured
Number (1,2)	Values	Values	Values
1	0.070	0.070	0.050
2	0.070	0.070	0.070
3	0.080	0.080	0.070
4	0.080	0.070	0.080
5	0.080	0.070	0.040
6	0.090	0.070	0.050
7	0.080	0.080	0.060
8	0.070	0.080	0.070
9	0.090	0.070	0.080
10	0.080	0.080	0.060
(1,2) Number Measurements (n)=	10	10	10
Standard Deviation (S)=	0.007	0.005	0.013
(4) Standard Conc. / S =	11	15	6
(5) Student's t ($\sigma = 0.99, n-1$) =	2.8214	2.8214	2.8214
(6) Lower Limit of Detection (LLD) =	0.02082	0.01457	0.03774
(6) Lower Limit of Quantitation (LLQ) =	0.10409	0.07285	0.18868
Lower Limit of Quantitation (μg) =	0.002 ug/filter	0.002 ug/filter	0.004 ug/filter



Savannah River

Nuclear Solutions, LLC

A Fluor Daniel Partnership

Method Reporting Limit for Ghostwipes™

Standard Conc. Used =

0.10

0.10

0.10

Units

ppb

ppb

ppb

Other Parameters/Information:

ghostwipe

ghostwipe

ghostwipes

Replicate

Measured

Measured

Measured

Number (1,2)

Values

Values

Values

1

0.090

0.090

0.100

2

0.100

0.090

0.100

3

0.100

0.090

0.100

4

0.100

0.100

0.080

5

0.100

0.100

0.100

6

0.090

0.130

0.100

7

0.110

0.090

0.100

8

0.100

0.100

0.090

9

0.100

0.090

10

0.110

0.100

^(1,2)Number Measurements (n)=

10

8

10

Standard Deviation (S)=

0.007

0.014

0.007

⁽⁴⁾Standard Conc. / S =

15

7

14

⁽⁵⁾Student's t ($\sigma = 0.99$, n-1) =

2.8214

2.9979

2.8214

⁽⁶⁾Lower Limit of Detection (LLD) =

0.01881

0.04066

0.01973

⁽⁶⁾Lower Limit of Quantitation (LLQ) =

0.09405

0.20329

0.09864

⁽⁷⁾Lower Limit of Quantitation (μg) =

0.009ug/wipe

0.020 ug/wipe

0.010 ug/wipe




Savannah River

Nuclear Solutions, LLC

A Fluor Daniel Partnership



Recovery of Ghostwipes™ + HP™ MCE BeO Standards



Recovery of Ghostwipes™ + HP™ MCE BeO Standards

Concentration	Expected Value µg	Result ppb	Result µg	% Recovery
BeO 0.01*	0.01	0.14	0.014	140.0
BeO 0.01*	0.01	0.11	0.011	110.0
BeO 0.2	0.2	2.01	0.201	100.5
BeO 0.2	0.2	1.92	0.192	96.0
BeO 0.5	0.5	4.88	0.488	97.6
BeO 0.5	0.5	4.51	0.451	90.2
BeO 1.0	1	9.19	0.919	91.9
BeO 1.0	1	8.97	0.897	89.7

*Concentration below established Method Reporting Limit.



Interference Testing, Fe/Ti/Pu in Ghostwipes™

Interference	Conc Interference	Conc Be	Be result	Recovery	Interference	Conc Interference	Conc Be	Be result	Recovery
Fe	10 mg	0.2	0.2	100	Ti	2 mg	0.2	0.21	105
	10 mg	0.2	0.16	125		2 mg	0.2	0.22	110
	10 mg	2	1.74	87		2 mg	2	1.72	86
	10 mg	2	1.76	88		2 mg	2	1.69	85
	10 mg	0.2	0.21	95	Pu	~5E6dpm	0.2	0.21	105
	10 mg	0.2	0.2	100		~5E6dpm	0.2	0.2	100
	10 mg	0.2	0.27	135		~5E6dpm	2	1.85	93
	10 mg	0.2	0.23	115	HP-CRMD	n/a	0.1	0.104	104
	10 mg	2	1.84	92	HP-CRME	n/a	0.2	0.201	101
	10 mg	2	1.8	90					
Blank	n/a	0	0.04	n/a					
Blank	n/a	0.01	0.04	n/a					
Blank	n/a	0.01	0.02	n/a					






Ghostwipe™ Lot Comparison:

Ghost Wipe Lot Comparison			
	Lot 1	Lot 2	Lot 3
	<0.005µg	<0.005µg	<0.005µg
	<0.005µg	<0.005µg	<0.005µg




Quality Control Chart of Beryllium Oxide by Optical Fluorescence


High Purity™ Certified Reference Material Beryllium Acetate Filters




HP Beryllium Oxide Filters @ 0.01 μ g Be



Ghostwipes™ Spiked with Beryllium Acetate



High Purity™ Filters Spiked with Beryllium Acetate





Summary

- **Savannah River Site currently uses Ghostwipes for all wet surface sampling. For Industrial Hygiene to convert to the optical fluorescence method, very little change would be required in the field.**
- **The addition of Optical Fluorescence as a Field of Testing allows the customer an option that should result in faster turnaround times and save money (labor dollars). The exact savings have not yet been calculated but are expected to be at least 50% of the current costs.**
- **The manual system is sufficient to handle SRNS current sample load, but if that changes, the option for an automated unit is available.**



Future Needs

- **There is no established standard for wet wipes used for beryllium surface sampling.**
- **Commercially available wipes are not required to adhere to a standard that ensures their suitability for beryllium measurement.**
- **Differences between lots in some brands of wipes is significant.**
- **There is no promise from any commercial provider that the wipes we use for beryllium sampling will not vary in material, preservative, or even background beryllium concentration, and this is a concern for ICP as well as Optical Fluorescence use.**



Further Information

Linda Youmans-McDonald
803-952-7885

linda.youmans-mcdonald@srs.gov

Mike Brisson
803-952-4400

mike.brisson@srs.gov

Anoop Agrawal
520- 321-7680

aagrawal@qwestoffice.net

John Cronin
520- 321-7680

jcronin@qwestoffice.net

Maureen Bernard
803-952-3500

maureen.bernard@srs.gov





Acknowledgements

It takes the cooperation of many to do this work. This slide attempts to recognize a few.

- Shirley Riley
- Erika Lee
- Marcia Gibson
- Matthew Hammett

Thanks! We could not have done it without you!!!

