




We Put Science To Work

Optical Fluorescence Analysis of Beryllium Wipes Contaminated with High Levels of Plutonium

Chad N. Durden, Linda Youmans-McDonald, and Michael J. Brisson

October 3, 2012

10:35AM



4th International Symposium on Beryllium Particulates and Their Detection

SRNS-L4600-2012-00174

Disclaimers

- **Mention of commercial products in this presentation does not imply endorsement by the author, SRNS, SRNL, SRS, or the U.S. Department of Energy (DOE).**
- **The speaker is an employee of a DOE contractor and is not a spokesperson for DOE itself.**

Summary


- **Introduction**
- **Fluorometer Glovebox Installation**
- **Method Testing**
- **Demonstration of Competency**
- **Interference Testing and Results (Ce and Pu)**
- **Conclusion**
- **Acknowledgements**

Introduction

- **Savannah River Site (SRS) radiological IH lab tasked to provide AIHA accredited Be results for GhostWipes™ contaminated with high amounts of plutonium (Pu)**
- **Up to 0.5 gram of Pu/wipe (samples from Pu processing glovebox)**
 - Customer assurance of <0.5 gram of Pu/wipe
- **Routine ICP-AES wipe samples undergo ion exchange (TEVA and Diphonix) for interference removal**
 - Modified NIOSH 7303 ICP, EPA 6010D
 - Glovebox modified ICP-AES is cost prohibitive, long lead time
- **Optical fluorometer and hotblock glovebox installation**
 - Existing method setup in radiohood
 - Comparable to ICP-MS sensitivity
 - More practical, cost effective, and less time to establish compared to a glovebox modified ICP-AES


Glovebox Fluorometer Installation

- **Introduce equipment into glovebox**
 - Fluorometer, hotblock, analytical balance, pipettes, thermometer, and associated cables (USB and power)
 - In-house addition of Amphenol connectors to USB cable for glovebox interfacing
 - Promega GloMax Multi Jr. fluorometer (FI-B) equipped with optical kit FM-B obtained from Berylliant



Glovebox Fluorometer Installation, continued

- **Remotely connect fluorometer console to PC**
- **Establish communication between PC and fluorometer**
- **Method development**
- **Resolve glovebox related, ergonomic issues**
 - Minimize unnecessary glovebox work, prepare as much as possible on benchtop
 - Utilize a click-pen with ink removed as stylus for fluorometer touchscreen
 - Readily available stylus' were too small for glovebox gloves.
 - Utilize long, plastic forceps for handling and manipulating cuvette caps




NIOSH 9110 and Other References

- **NIOSH 9110 – Beryllium in Surface Wipes by Field-portable Fluorometry**
 - Alternative to NIOSH 7300 (hotblock digestion and ICP-AES analysis)
 - **Purpose: accurately measure beryllium using a basic solution to remove metal contaminants and a dye specific to the beryllium atom.**
 - Hydroxybenzoquinoline sulfonate (HBQS) specific to beryllium
 - Ethylenediaminetetraacetic acid (EDTA) and high alkalinity remove interferences. e.g., Fe and Ti
- **DOE action limit 0.2 $\mu\text{g}/100 \text{ cm}^2$**
- **Implement as much as possible from existing publications and SRS radiohood experience**
 - Ashley, K.; Agrawal, A.; Cronin, J.; Tonazzi, J.; McClesky, T. M.; Burrell, A. K.; Ehler, D. S. Ultra-trace determination of beryllium in occupational hygiene samples by ammonium bifluoride extraction and fluorescence detection using hydroxybenzoquinoline sulfonate. **Analytica Chimica Acta** 584. 2007, 281-286.
 - ASTM D7202 - Standard Test Method for Determination of Beryllium in the Workplace Using Field-Based Extraction and Fluorescence Detection

Method Testing

- **Interference testing conducted with and without blank sampling media (GhostWipe™)**
- **3 wt% ammonium bifluoride (ABF) selected for glovebox wipes**
 - Assured dissolution of refractory BeO in highly soiled wipe samples
 - 1 wt% ABF routinely used for wipes and filters in the radiohood
- **Establish working calibration prepared in 3 wt% ABF**
 - 0.00, 0.25, 0.50, 2.0, and 10.0 ppb at instrument
 - Linear calibration range: $R^2 \geq 0.998$
- **BeO CRM filter extractions in 3 wt% ABF**
 - Obtained from High Purity Standards
 - Prepared from NIST SRM 1877 (high-fired BeO powder)

| Calibration Data | |
|------------------------------|-------------------|
| Standard Concentration (ppb) | Instrument Counts |
| 0.00 | 108.80 |
| 0.25 | 307.20 |
| 0.50 | 485.80 |
| 2.00 | 1590.40 |
| 10.00 | 7533.30 |



Method Testing, continued

- **Wipes submerged in 20 mL of 3 wt% ABF**
- **Heated at 85-90°C for 1 hour**
- **Extractions allowed to cool and settle overnight**
 - Care was taken not to disturb extracted solutions during cooling and aliquoting
- **1 mL of extract was filtered through a 0.2 micron filter**
- **100 µL of filtered extract was added to 1.9 mL of detection solution**
 - Aliquots were taken near the surface of solution
 - Total 400x dilution (20x extraction + 20x detection)
- **Preparations were placed in darkness ~2 hours**
- **Solutions were measured by fluorometry**
 - pH confirmed >12

Demonstration of Competency


- **Demonstration of Competency (DoC) tests are performed in-house in place of Proficiency Analytical Test Samples (PATS)**
 - AIHA BePATS are analyzed by ICP-AES
 - DoC consists of CRM filters blind to the analyst
- **Six BeO CRM filters extracted in 3 wt% ABF**
 - Three 0.2 µg BeO/filter CRMs
 - Three 0.5 µg BeO/filter CRMs
- **Average BeO recovery = 92%**
 - AIHA acceptable recovery range 75 – 125%

| | Reference (µg/filter) | Measured (µg/filter) | % recovery |
|---------|-----------------------|----------------------|------------|
| BeO 0.2 | 0.2 | 0.172 | 86.0% |
| BeO 0.2 | 0.2 | 0.192 | 96.0% |
| BeO 0.2 | 0.2 | 0.160 | 80.0% |
| | | | |
| BeO 0.5 | 0.5 | 0.472 | 94.4% |
| BeO 0.5 | 0.5 | 0.520 | 104.0% |
| BeO 0.5 | 0.5 | 0.460 | 92.0% |

Plutonium and Optical Fluorescence

- **Limited data available describing the effect of Pu on Be by fluorometry**

- Fluoride necessary to ensure complete dissolution of BeO



- Fluoride could be complexed with Pu, not allowing complete dissolution of BeO
 - A small amount of HF (0.05-0.1M HF/8-12M HNO₃) is desirable for PuO₂ dissolutions
- Complicated optical spectra of actinides justified PuO₂ interference testing



Interference Testing

- **CeO₂ used in initial testing**
 - Ce is safer, cheaper, easier to obtain, and easier to dispose than Pu
 - Obtained from Aldrich, >99.9% metals basis, <5 micron powder
 - CeO₂ was utilized instead of Ce metal to minimize H₂ gas generation
 - Ce was used in cold runs at SRS before Pu processing Canyons were placed online
- **PuO₂ used to confirm CeO₂ testing**
 - NBL CRM No. 122 PuO₂ in powder form
 - Limited to 3 grams of CRM No. 122
 - <1 µg/g beryllium impurity
 - Concentrated Pu solutions were not considered due to low pH

New Brunswick Laboratory Certified Reference Materials Certificate of Analysis

CRM No. 122
Plutonium Oxide - PuO₂
In Powder Form
(Plutonium Assay and Isotopic Standard)
(In cooperation with the University of California
Los Alamos National Laboratory, Los Alamos, New Mexico)

| | |
|------------------------------|--|
| Plutonium | 87.790 ± 0.039 Wt. %* (877.90 ± 0.39 g·kg ⁻¹) |
| Plutonium-238 | 0.0521 ± 0.0011 At. %* |
| Plutonium-239 | 87.305 ± 0.004 At. %* |
| Plutonium-240 | 11.539 ± 0.004 At. %* |
| Plutonium-241 | 0.9248 ± 0.0011 At. %* |
| Plutonium-242 | 0.1790 ± 0.0013 At. %* |
| Relative Atomic Weight | 239.191* |

Interference Testing – Cerium Oxide

- **Performed with and without blank sampling media (GhostWipe™)**
- **3 wt% ABF solvent**
 - CeO₂ used as PuO₂ surrogate
- **CeO₂ material tested at:**
 - 0.1 gram CeO₂
 - 0.25 gram CeO₂
 - 0.5 gram CeO₂
- **Undissolved CeO₂ material was observed (pale yellow) following heat cycle**
- **ABF solution was colorless following heat cycle**

| Sample ID | Instrument Result (ppb) | |
|---------------------------|-------------------------|----------------------|
| | with blank wipe | without blank wipe |
| Prep Blank (PB) | 0.03 | |
| Indep. Cal. Verif. (ICV) | | 2.00 |
| BeO 0.2 µg Be | | 0.47 (0.19 µg Be) |
| BeO 0.5 µg Be | | 1.33 (0.53 µg Be) |
| BeO 0.2 µg Be | 0.43 (0.17 µg Be) | |
| BeO 0.5 µg Be | 1.18 (0.47 µg Be) | |
| LCS 0.25 ppb Be, no Ce | 0.24 | 0.27 |
| LCSD 0.25 ppb Be, no Ce | 0.26 | 0.42 |
| LCS 0.5 ppb Be, no Ce | 0.47 | 0.59 |
| LCSD 0.5 ppb Be, no Ce | 0.46 | 0.52 |
| LCS 2.0 ppb Be, no Ce | 1.78 | 2.04 |
| LCSD 2.0 ppb Be, no Ce | 1.78 | 2.11 |
| LCS 0.25 ppb Be, 0.1g Ce | 0.24 | 0.24 |
| LCSD 0.25 ppb Be, 0.1g Ce | 0.26 | 0.35 |
| LCS 0.5 ppb Be, 0.25g Ce | 0.45 | 0.53 |
| LCSD 0.5 ppb Be, 0.25g Ce | 0.45 | 0.53 |
| LCS 2.0 ppb Be, 0.5g Ce | 1.72 | 2.05 |
| LCSD 2.0 ppb Be, 0.5g Ce | 1.74 | 2.07 |

Interference Testing – Plutonium Oxide

- Performed with and without blank sampling media (GhostWipe™)
- 3 wt% ABF solvent
 - PuO₂ powder used to confirm CeO₂ testing
- PuO₂ CRM material tested at:
 - 0.1 gram PuO₂
 - 0.25 gram PuO₂
 - 0.5 gram PuO₂
- Undissolved PuO₂ material was observed (black) following heat cycle
- ABF solution was colorless following heat cycle

| Sample ID | Instrument Result (ppb) | |
|------------------------|-------------------------|--------------------|
| | with blank wipe | without blank wipe |
| Prep Blank (PB) | 0.03 | |
| Indep. Cal. Ver. (ICV) | | 1.99 |
| | | |
| BeO 0.2 µg Be | 0.48 (0.19 µg Be) | |
| BeO 0.5 µg Be | 1.30 (0.52 µg Be) | |
| | | |
| 0.5g Pu/no Be | 0.03 | 0.04 |
| 0.25g Pu/no Be | 0.05 | 0.04 |
| 0.25g Pu/no Be | 0.03 | 0.02 |
| 0.1g Pu/no Be | 0.11 | 0.04 |
| 0.1g Pu/no Be | 0.03 | 0.02 |

Conclusion

- **Optical fluorescence method successfully setup in glovebox**
- **CeO₂ and PuO₂ interference testing confirmed no significant impact to beryllium measurements at tested levels**
- **Method utilized on “real” samples**
- **AIHA accredited Be results were reported to the customer**
- **Lower report limit = 0.1 µg Be/wipe**
 - 0.25 ppb at instrument
 - Total 400x dilution (20x extraction + 20x detection)
 - Additional testing necessary to obtain lower report limit

Acknowledgements

- **Sincere gratitude is extended to lab technician Shirley W. Riley for her patience, dedication, and commitment to this project.**
- **Gratitude is extended to Erika M. Lee.**
- **Gratitude is also extended to Kevin Ashley, Berylliant, and LANL.**



Questions?