

The Beryllium Lymphocyte Proliferation Test (BeLPT)

Lisa Barker
National Jewish Health



Recognition of Beryllium Health Effects

1940s: Reports of adverse health effects associated with beryllium exposure.

- Acute disease with high exposures.
- Chronic disease present.

1950s: Be exposure standard ($2 \mu\text{g}/\text{m}^3$) implemented but cases of CBD continue to be reported

Recognition of Beryllium Health Effects

1980s: New test (BeLPT) identifies beryllium sensitization, a precursor to CBD.

2000: Department of Energy mandates CBD protection programs

2001: EEOICPA: Compensates employees injured during employment in the nuclear weapons development effort

Extent Of Beryllium Use

- NIOSH estimates
 - At risk: 134,000 in the US
(J Occup Environ Hyg 2004 1:648)
 - Probably conservative- Doesn't include “downstream users”
- Study of Construction Trades Workers at 3 DOE sites
 - Estimate 230,000 exposed (Welch 2004)





Importance of an Accurate Diagnosis

- Intervene early
 - If CW, remove/reduce exposure
 - Reduce exposures to workers in same job
- If CBD
 - Follow closely for progression of symptoms
 - Advise to maintain health (exercise, preventive vaccines)
 - Treat to preserve lung function

Beryllium Lymphocyte Proliferation Test (BeLPT)



- Measures blood cells response to beryllium
- Identifies beryllium sensitization




- Further clinical testing needed to diagnose CBD

BeLPT Characteristics

- Detects a beryllium-specific cell mediated immune response Identify CBD at early stages
- Differentiates CBD from lung disease of other etiology (sarcoidosis, COPD, HP)
- Increased sensitivity and specificity to diagnose BeS/CBD compared with other tests
 - Occupational History
 - CXR or CT scan
 - Spirometry




The BeLPT is Specific to Be - Related Health Effects



How is the Test Performed

- Culture peripheral blood or BAL cells with and without beryllium
- Measure incorporation of ^3H -Td for proliferation
- Stimulation index or SI : ratio of Be-stimulated to unstimulated cpm
- Use positive controls to ensure cells are able to proliferate



Lymphocyte Proliferation Test: Peripheral Blood

NORMAL LPT

Patient Name:

ID Number:

OS#:

Date of Test:

Accession#:

Referring Physician:

Technician:

Results

Mean Stimulation Index

Mitogens:

Day 3

Day 4

Day 6

Phytohemagglutinin

431.3

Antigens:

Candida

107.8

Beryllium Sulfate:

1X10⁻⁴ M

1.1

0.6

1X10⁻⁵ M

1.5

0.8

1X10⁻⁶ M

1.1

0.6

Normal response to mitogen.

Normal response to antigen. (Normals: >3.0)

Normal lymphocyte proliferation to beryllium sulfate.

Note: A normal blood result does not exclude the diagnosis of Chronic Beryllium Disease. If the clinical setting is compatible, a beryllium lymphocyte proliferation test using bronchoalveolar lavage lymphocytes is recommended.

A normal result is all values equal or below the cutoff value of 2.5.

Raw Data for Normal Blood LPT

| PHA | | CANDIDA | | DAY 4 | | | | DAY 6 | | | | |
|-------|---------|---------|---------|-------|-------|-------|-------|--------|-------|-------|-------|------|
| | | | | C+M | | | | C+M | | | | |
| | | | | _342 | | | | _393 | | | | |
| | | | | _215 | | | | _560 | | | | |
| | | | | _243 | | | | _351 | | | | |
| | | | | _197 | | | | _742 | | | | |
| | | | | _221 | | | | _455 | | | | |
| | | | | _336 | | | | _545 | | | | |
| | | | | _190 | | | | _603 | | | | |
| | | | | _159 | BER-4 | BER-5 | BER-6 | _884 | BER-4 | BER-5 | BER-6 | |
| C+M | PHA | C+M | CAN | 675* | _229 | _452 | _217 | 516 | _402 | _514 | _330 | |
| _179 | _88819 | _452 | _93302* | 403* | _283 | _399 | _256 | 801 | _312 | _549 | _433 | |
| _288 | _97529 | _423 | _49888 | 430* | _237 | _301 | _364 | 799 | _296 | _382 | _379 | |
| _202 | _82100 | _196* | _46436 | 278 | _322 | _256 | _276 | 919* | _154* | _386 | _310 | |
| _162 | _89930 | _629 | _65880 | 242.3 | 267.8 | 352.0 | 278.3 | 604.5 | 336.7 | 457.8 | 363.0 | MEAN |
| 207.8 | 89594.5 | 501.3 | 54068.0 | 64.1 | 43.3 | 89.5 | 62.2 | 178.4 | 57.1 | 86.4 | 54.9 | SD |
| 56.0 | 6320.3 | 111.5 | 10374.1 | 0.26 | 0.16 | 0.25 | 0.22 | 0.30 | 0.17 | 0.19 | 0.15 | CV |
| 0.27 | 0.07 | 0.22 | 0.19 | | 1.1 | 1.5 | 1.1 | | 0.6 | 0.8 | 0.6 | SI |
| | 431.3 | | 107.8 | | | | | | | | | |
| | | | | | | | | | | | | |
| Norm: | | Norm: | | | | | | Normal | | | | |

Calculate Stimulation Index (D4B4) $\frac{267.8}{242.3} = 1.1$

Lymphocyte Proliferation Test: Peripheral Blood

ABNORMAL LPT

Patient Name:

ID Number:

OS#:

Date of Test:

Accession#:

Referring Physician:

Technician:

Results

Mean Stimulation Index

Mitogens:

Day 3

Day 4

Day 6

Phytohemagglutinin 125.5

Antigens:

Candida

172.6

Beryllium Sulfate:

1X10⁻⁴ M

8.0

39.0

1X10⁻⁵ M

5.0

30.0

1X10⁻⁶ M

1.9

2.7

Normal response to mitogen.

Normal response to antigen. (Normals: >3.0)

Abnormal lymphocyte proliferation to beryllium sulfate.

Note: An abnormal result is two or more values above the cutoff value of 2.5.

Raw Data for Abnormal LPT

| PHA | | CANDIDA | | DAY 4 | | | | DAY 6 | | | | |
|-------|---------|---------|----------|--------|--------|--------|--------|---------|---------|---------|--------|------|
| | | | | C+M | | | | C+M | | | | |
| | | | | _679 | | | | _811 | | | | |
| | | | | _1246 | | | | _739 | | | | |
| | | | | _919 | | | | _713 | | | | |
| | | | | _1012 | | | | _523 | | | | |
| | | | | _1136 | | | | _479 | | | | |
| | | | | _370* | | | | _837 | | | | |
| | | | | _529 | | | | _750 | | | | |
| C+M | PHA | C+M | CAN | _1918* | BER-4 | BER-5 | BER-6 | _274 | BER-4 | BER-5 | BER-6 | |
| _580 | _81214 | _1394* | _121866 | _75* | _6349 | _8521* | _1460 | _429 | _40508* | _18750 | _3327* | |
| _648 | _64481 | _892 | _103238 | _851 | _8206 | _4178 | _1838 | _681 | _21183 | _5471* | _2368 | |
| _368 | _66263 | _588 | _103039 | _777 | _7032 | _5936 | _1809 | _1067* | _30678 | _15870 | _1511 | |
| _509 | _52115 | _403 | _105203 | _1662* | _7132 | _3264 | _1833 | _1564* | _21077 | _21594 | _1259 | |
| 526.3 | 66018.3 | 627.7 | 108336.5 | 893.6 | 7179.8 | 4459.3 | 1735.0 | 623.6 | 24312.7 | 18738.0 | 1712.7 | MEAN |
| 119.8 | 11925.2 | 246.9 | 9072.4 | 236.6 | 767.6 | 1358.0 | 183.8 | 186.3 | 5512.8 | 2862.0 | 581.4 | SD |
| 0.23 | 0.18 | 0.39 | 0.08 | 0.26 | 0.11 | 0.30 | 0.11 | 0.30 | 0.23 | 0.15 | 0.34 | CV |
| | 125.5 | | 172.6 | | 8.0 | 5.0 | 1.9 | | 39.0 | 30.0 | 2.7 | SI |
| | | *HI CV* | | | | | | | | | | |
| Norm: | | Norm: | Keep | | | | | Abnorm: | | | | |

Calculate Stimulation Index (D6B4) $24312.7 / 623.6 = 39.0$

Lymphocyte Proliferation Test: Peripheral Blood

UNINTERPRETABLE

Patient Name:

ID Number:

OS#:

Date of Test:

Accession#:

Referring Physician:

Technician:

Results

Mean Stimulation Index

Mitogens:

Day 3

Day 4

Day 6

Phytohemagglutinin 36.0

Antigens:

Candida

4.0

Beryllium Sulfate:

1X10⁻⁴ M

1X10⁻⁵ M

1X10⁻⁶ M

Normal response to mitogen.

Normal response to antigen. (Normals: >3.0)

Uninterpretable lymphocyte transformation to beryllium sulfate.

Why is a Test Uninterpretable

- Cells “over” respond in all conditions, possibly masking a true lymphocyte proliferation in Be stimulated condition(s)
- Positive controls not positive
- Lab error

Specificity and Sensitivity of the BeLPT

Maier AOEH, 2001

| Study | Individuals Enrolled, n | BeS/CBD, n (%) | Estimated Sensitivity | Estimated Specificity | Positive Predictive Value |
|---|-------------------------|----------------|-----------------------|-----------------------|---------------------------|
| Mroz et al, JACI 1991 | 35 | 17 (49%) | 94% | 100% | |
| Stokes et al, JOEM 1991 | 57 | 27 (47%) | 38% | 97% | |
| Kreiss et al, JOEM 1993 ⁺ * | 505 | 9 (2%) | 89% | NA | 100% |
| Kreiss et al, ARRD 1993 ⁺ * | 895 | 18 (2%) | 100% | NA | 100% |
| Stange et al, Environ Health Perspect 1996 ⁺ * | 4268 | 101 (2%) | 97% | NA | |


Predictive Value of Clinical Tests in the Diagnosis of CBD

Kreiss et al JOEM 1989

| Test | Positive Predictive Value | Negative Predictive Value |
|-----------------------|------------------------------|------------------------------|
| Symptoms | 1.6-10.5% | 98.3-99% |
| FEV ₁ | 2.9% | 98.3% |
| FVC | 3.3% | 98.2% |
| FEV ₁ /FVC | 3.4% | 98.5% |
| Chest Xray | 12.9% | 98.9% |
| LPT | 100% | 99.8% |

Alternative Screening Tests for CBD

- Chest X-ray
 - PV+ = 12.9% (Kreiss et al, 1993)
 - PV+ = 0.8% (Stange et al, 1996)
- Spirometry
 - PV+ of FEV_1 = 2.9% (Kreiss et al, 1993)
 - PV+ of FVC = 3.3% (Kreiss et al, 1993)
- Symptoms
 - PV+ of COUGH = 3.2% (Kreiss et al, 1993)
 - PV+ of DYSPNEA = 5.0% (Kreiss et al, 1993)



Interlaboratory Agreement

- Previous report 85-96% between labs, 30-60% for abnormal tests (Stange, Env Health Perspec 1996)
- Prelim evaluation split test from 2 labs (n=825)
 - 716/825 agreed (87%)
 - 703/716 negative (98% of agreed, 85% of total)
 - Total positive: 44 (5%)

BeLPT compared to other medical screening tests

| | Test | Sensitivity | Specificity | Follow-up test |
|---------------------|--------------------|-------------|-------------|-------------------------------|
| Breast CA | Mammography | 75 – 88% | 83 – 99% | Biopsy |
| Colon CA | Fecal Occult Blood | 26 – 92% | 90 – 99% | Colonoscopy |
| Neural Tube Defects | AFP | 56 – 91% | 5 – 10% | Amniocentesis |
| Prostate CA | PSA | 29 – 80% | 75% | Biopsy |
| BeS/CBD | BeLPT* | 74 – 91% | 96 – 99% | Bronchoalveolar lavage/biopsy |

Considerations in using the BeLPT

- Borderline/uninterpretable
- False negative
- False positives
- Lack of concordance between some labs
- At least 7 days for results
- Does not discriminate between BeS and CBD

Efforts to Standardize Testing

- CABST
- ICPT
- DOE Specification, April 2001 (update expected)
 - Minor differences remain:
 - Days of culture
 - Cut-off value determination
 - Serum lots



What if There are Problems with the BeLPT?

- Sometimes the cells fail to respond to beryllium. A negative BeLPT does not exclude the diagnosis of CBD
 - More common in smokers
- Beryllium skin patch testing can be used to demonstrate an immune response to beryllium




Benefits of the BeLPT

- Redefined and broadened our understanding of beryllium related health effects:
- Can detect health effects at an early stage
- Higher sensitivity and specificity than other screening tests
- Can be used in screening and surveillance mode to help evaluate processes/tasks

Applications of the BeLPT

- Diagnosis of beryllium related health effects
 - Beryllium Sensitization
 - CBD
- Workplace screening and surveillance
 - Prevalence of health effects
 - Work practices and job titles with increased risk



Summary

- BeLPT provides non-invasive test to screen for Be-health effects
 - Specific and sensitive
 - Used in evaluation of process/exposure risks
- Studies needed of other tests
 - Discriminate BeS and CBD
 - Enhance performance of BeLPT
 - Increased specificity/sensitivity

Questions ?

